

# Cool Drape Too

MAN & MACHINE

Maintaining a hygienic workspace has never been more accessible with our premium protective keyboard cover. Designed to fit any keyboard, it can be placed on top or wrapped around for full coverage.

Made from ultra-thin, latex-free silicone, Cool Drape Too is transparent, ensuring clear key visibility while providing a reliable barrier against germs and contaminants. Cleaning is effortless - simply rinse it under a tap, immerse it in water, expose it to UV disinfection, or even sterilise it in an autoclave. Durable, reusable, and practical, Cool Drape Too is the perfect choice for any setting where hygiene is a priority.



## CLEAN HYGIENICALLY WASHABLE

To fit your cleaning protocol, the drape may both be dipped, rinsed, sprayed, or wiped with soap and water as well as approved cleaning fluids.



## USE YOUR STANDARD DISINFECTANT DISINFECTABLE

The drape resists most common surface disinfectants, allowing you not only to clean it, but also disinfect it. Treat it just like any other high frequency touch surface in your facility.



## ULTIMATE CLEAN STERILISABLE

Designed for the highest level of sanitation, this drape can be repeatedly sterilised in an autoclave or washer-disinfector at temperatures up to 134°C.



## DURABLE REUSABLE

Built for longevity, the drape withstands repeated sanitation and can be stored or reapplied for ongoing use.



## TOP PERFORMANCE KEEP CONTROL

A slip-resistant design ensures smooth and reliable access to keys and the touchpad for secure data entry precision.



## KEY CONTROL OPTIMAL TRANSPARENCY

The ultra-thin silicone sheet is designed for exceptional transparency while offering a soft, smooth, and responsive feel for a comfortable typing experience.

# Cool Drape Too

## Cleaning Specifications

### CLEANING FLUIDS ALLOWED



Soap and water.



Approved surface disinfectants, such as those based on alcohol, chlorine  $\leq 10\%$ , formaldehyde  $\leq 37\%$ , glutaraldehyde, phenol or quaternary ammonium  $\leq 2\%$ , or oxidising  $\leq 3$ .

### CLEANING & STERILISATION METHODS ALLOWED



You can wipe the product with a soft cloth that is non-abrasive like a microfiber cloth, a disinfecting wipe, a soft brush similar to a paint brush in bristle stiffness or a sponge, but not a scourer of any sort.



You can spray the product.



You can rinse the product.



You can immerse the product.



You can expose the product to Hydrogen Peroxide Vapour (HPV).



You can sterilise the product in an autoclave or thermo disinfectant (up to 134 °C, tested 200 cycles).



You can expose the product to UV light disinfection.

### FULL LIST OF APPROVED SURFACE DISINFECTANTS



### PRODUCT MANUAL



## Product Specifications

**PRODUCT DESCRIPTION** Protective covers to use with any keyboard, pack of 4

**MATERIAL COMPOSITION** Latex-free silicone

**WARRANTY** No manufacturer warranty applies

**COMPLIANCE** CE

**COLOUR** Transparent

**SIZE** 45 x 25 cm

**WEIGHT** net. 50 g, 200 g per 4-pack

**STORAGE TEMPERATURE** -20 °C to 70 °C

**OPERATING TEMPERATURE** 0 °C to 70 °C; Sterilisable at 134 °C

**ORIGIN** Made to our specifications in Taiwan

